

www.gmmspa.com







STATE-OF-THE-ART DESIGN AND TECHNOLOGY

Adjustable-height examination table with tabletop floating in the four directions: easy and fast patient positioning and maximum patient throughput.





The compact and solid structure of the ceiling stand provides flexible and smooth movements.

- Manual or automatic collimation
- Auto-tracking
- Wide touch screen ensuring the control of the complete system



An advanced technology provides the automatic alignment of the X-ray source with the detector movement so that the operator can easily execute AP projections and oblique incidences even varying the SID.



The wall stand can be vertically adjusted and tilted fot the best execution of examinations of the skull, thorax, column and limbs.



The detector 90°/20° tilt from the vertical position facilitates patient positioning for special incidences of the skull and joints. The detector 90° rotation from the horizontal position ensures easy and quick examinations of patients carried on stretchers.



MAXIMUM OPERATIONAL RELIABILITY









